


# Scans



For:  
Temporary prostheses on implants  
Provisional Denture




## Temporary prostheses on implants



# Case n°1

Preliminary - Provisional Scans  
One Shot IQ with soft tissues



extra 2: reverse scan  
with provisional  
(optional for lower jaw)

extra 1: one shot with  
soft tissues




## Provisional Denture



# Case n°1


preliminary - with relined denture



Mesh extra 1: one shot IQ with soft tissues

# Case n°2


Scan Flag IQ Scans



Recommended for Scan flags IQ

# Tutorial

Scan of Scan Flag IQ



[Play Video](#)